



## TECHSPEC® VARIABLE BEAM EXPANDER INSTRUCTIONS

## 1X-3X, 266nm, #87-559



## Variable Beam Expander Magnification Settings 1X-3X, 266nm (#87-559)



Beam expander magnification and collimation is set by adjusting both the Alpha and Numeric rings using the magnification settings chart as a guideline. Nominal magnification is indicated along the x-axis, and the settings of the Alpha and Numeric rings are indicated along the y-axes. The **blue curve** follows the rotation of the Alpha ring and the **red curve** follows the rotation of the Numeric ring. **Note:** These are nominal settings. Some adjustment for your individual application may be necessary for optimal performance.

To set the beam expander to a specific magnification:


- Find the desired nominal magnification along the x-axis and trace a vertical line up through the blue and red curves.
- 2. The intersection of that vertical line with the **blue curve** corresponds to a letter on the left y-axis. Rotate the alpha ring to that location. Hash marks are provided between each letter on the ring to help dial in the required magnification.
- The intersection of that vertical line with the **red curve** corresponds to a number on the right y-axis. Rotate the numeric ring to that location. Hash marks are provided between each number on the ring to help dial in the required magnification.

**For Example -** To set a magnification of 2.6X, turn the alpha ring to G and the numeric ring to 9.

part #84132110

www.edmundoptics.com/vari-beam

For Stock No. 87-559



www.edmundoptics.com/vari-beam

part #84132110

For Stock No. 87-559