

PHOTONICS WORKSHOP

LIMITLESS OPTIONS WITH EO STOCK PARTS

NOT SURE IF STOCK LASER OPTICS ARE THE RIGHT FIT FOR YOUR APPLICATION?

This demo shows how easy it is to integrate stock laser optics into R&D laser systems including Michelson Interferometers, beam combiners, and systems designed to characterize beam profiles. The options are limitless for R&D, laboratory, and prototyping applications. Looking for stock optics, optomechanics, and lasers? Think Edmund Optics®.

FEATURED STOCK PRODUCTS

TECHSPEC® Laser Optics

- 1750 Unique Laser Optics In Stock and Available for Immediate Delivery
- Lenses, Filters, Mirrors, and More

TECHSPEC® Optical Cage System

- Limitless Customization and Adaptability
- Ideal for Prototyping

TECHSPEC® Mechanics

- Easily Interchange TECHSPEC® Components
- Modular Mounting, Translation Stages, Slides, and More

Contact Us to Design Your Lab!

GERMANY: +49 (0) 721 6273730 UK: +44 (0) 1904 788600

FRANCE: +33 (0) 820 20 75 55 ITALY: +39 800 875 211

EO Edmund
optics | worldwide

www.edmundoptics.eu

READY TO BUILD YOUR OWN?



FEATURED STOCK PRODUCTS:

TECHSPEC® LASER OPTICS

www.edmundoptics.com/laser-optics

Description	Stock Number
TECHSPEC® 427nm Dichroic Laser Beam Combiner, 25mm Diameter	#86-389
TECHSPEC® 1/40λ Mirror, 25mm Dia, Enhanced Aluminum Coating	#64-015
TECHSPEC® 5X Beam Expander 532nm	#59-126
TECHSPEC® 532nm Laser Non-Polarizing Plate Beamsplitter 50mm x 50mm	#48-220
Microlens Array 10 x 10mm, 500µm Pitch, 1.2° Divergence	#64-479
TECHSPEC® Polarizing Cube Beamsplitter 25mm VIS	#49-002
TECHSPEC® C-Mounted Non-Polarizing Cube Beamsplitter, VIS Coating	#49-682
Axicon 140° Apex Angle, 25.4mm Dia, Uncoated	#83-781
TECHSPEC® Circular Variable Filter 50mm Diameter Mounted 0.04-3.0 OD	#54-536

TECHSPEC® CAGE SYSTEM

www.edmundoptics.com/cage-system

Description	Stock Number
TECHSPEC® Cage Sphere with (5) 30mm Ports and (1) 43mm Port	#85-624
TECHSPEC® 25mm Standard Cage Plate with 1/4-20 for Post Mounting	#85-633
TECHSPEC® 25mm Basic Cage Plate	#88-196
TECHSPEC® 25.4mm Kinematic Mount	#85-695
TECHSPEC® 25/25.4mm Cube Beamsplitter Mount	#85-688

TECHSPEC® OPTOMECHANICS

www.edmundoptics.com/optomechanics

Description	Stock Number
TECHSPEC® 1" x 1" Kinematic Mount, 2 Screws	#58-857
TECHSPEC® Post Holder, 2" Length, 1/4-20 Thread	#58-978
TECHSPEC® 40mm, Side Drive, Solid Top, Translation Stage	#66-406