

[See all 7 Products in Family](#)

2X Mitutoyo BD Plan Apo Objective

See More by [Mitutoyo](#)



Mitutoyo Brightfield and Darkfield (BD) Infinity Corrected Objectives

Stock **#18-402** **3 In Stock**

⊖ 1 ⊕ €1.290⁰⁰

ADD TO CART

| Volume Pricing | |
|----------------|-------------------------------|
| Qty 1+ | €1.290,00 each |
| Need More? | Request Quote |

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General

378-831-13 **Model Number:**

Compatible Tube Lens Focal Length (mm):
Focal Length: 200mm

Microscope Objective **Type:**

Infinity Corrected **Style:**

Manufacturer: Mitutoyo

Note: #18-653 or #18-652 are required for darkfield operation.

Physical & Mechanical Properties

Length excluding Threads (mm): 61.00

Maximum Diameter (mm): 44

Weight (g): 340

Optical Properties

Compatible Cover Glass Thickness (mm): N/A

Focal Length FL (mm): 100.00

Magnification: 2X

Numerical Aperture NA: 0.06

Resolving Power (μm): 5

Depth of Field (μm): 91.00

Field of View, 24 Diameter Field Eyepiece (mm): 12

Working Distance (mm): 34

Wavelength Range (nm): 435 - 655

Parfocal Length (mm): 95

Immersion Liquid: N/A

Sensor

Maximum Sensor Format: 2/3"

Threading & Mounting

Mounting Threads: M40 x 36 TPI

Environmental & Durability Factors

Operating Temperature ($^{\circ}\text{C}$): +5 to +40

Regulatory Compliance

RoHS 2015: Exempt

Certificate of Conformance: View

REACH 241: Contains SVHC(s)

Product Details

- Ideal for Brightfield and Darkfield Microscopy Applications
- Long Working Distances and High Numerical Apertures
- Plan Apochromat Design

Mitutoyo Brightfield and Darkfield (BD) Infinity Corrected Objectives offer long working distances with high numerical apertures. The high-quality plan apochromatic design allows these objectives to maintain high numerical apertures and contrast while correcting chromatic aberration in the visible wavelength region. These objectives feature M40 x 36 TPI mounting threads. Mitutoyo Brightfield and Darkfield (BD) Infinity Corrected Objectives are optimized for both brightfield and darkfield illumination for applications including surface inspection and advanced life science microscopy.

Technical Information

Mitutoyo BD Plan Apo Objectives
Typical Transmission



| Stock No. | BD Plan Apo | A | B | C | D | E | F | W.D. |
|-----------|----------------|------|------|------|------|------|------|------|
| #18-402 | 2X | 61.0 | 4.0 | 40.0 | 37.0 | 42.0 | 44.0 | 34.0 |
| #18-403 | 5X | 61.0 | 4.0 | 40.0 | 37.0 | 42.0 | 44.0 | 34.0 |
| #18-404 | 7.5X | 61.0 | 4.0 | 40.0 | 37.0 | 42.0 | 44.0 | 34.0 |
| #18-405 | 10X | 61.0 | 4.0 | 40.0 | 37.0 | 42.0 | 44.0 | 34.0 |
| #18-406 | 20X | 75.0 | 7.0 | 40.0 | 32.0 | 42.0 | 44.0 | 20.0 |
| #18-407 | 50X | 82.0 | 8.5 | 40.0 | 32.0 | 42.0 | 44.0 | 13.0 |
| #18-408 | 100X | 89.0 | 10.0 | 40.0 | 31.5 | 42.0 | 44.0 | 6.0 |
| Stock No. | BD Plan Apo SL | A | B | C | D | E | F | W.D. |
| #18-409 | 20XSL | 64.5 | 7.5 | 40.0 | 31.0 | 42.0 | 44.0 | 30.5 |
| #18-410 | 50XSL | 75.0 | 7.0 | 40.0 | 32.0 | 42.0 | 44.0 | 20.0 |
| Stock No. | BD Plan APO HR | A | B | C | D | E | F | W.D. |
| #18-411 | 50XHR | 89.9 | 4.8 | 40.0 | 32.0 | 42.0 | 44.0 | 5.2 |
| #18-412 | 100XHR | 93.7 | 4.7 | 40.0 | 32.0 | 42.0 | 44.0 | 1.3 |



Units: mm