

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

0.1 - 4.0 OD, 75 X 75mm, Kodak Wratten 2 Neutral Density Filter Kit (13 Filters)



Stock #15-115 [CONTACT US](#)

- 1 + €2.395,00

ADD TO CART

Qty 1+

€2.395,00

Volume Pricing

[Request Quote](#)

 Prices shown are exclusive of VAT/local taxes

Get Product Downloads



Dimensions (mm):
75 x 75 (Nominal)

Filters Included :
Optical Densities: [0.1](#), [0.2](#), [0.3](#), [0.4](#), [0.5](#), [0.6](#), [0.7](#), [0.8](#),
[0.9](#), [1.0](#), [2.0](#), [3.0](#), [4.0](#)

Length (mm):
75.00

Type:
Neutral Density Filter

Width (mm):
75.00

Number of Filters:
13

Regulatory Compliance

RoHS:
Compliant

Certificate of Conformance:
[View](#)