

In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union

SUBMIT

TECHSPEC® 138mm Telecentric Backlight Illuminator



138mm Telecentric Backlight Illuminator, #35-427

Stock #35-427 3-4 DAYS



1



€2.205,00

ADD TO CART

Qty 1+

€2.205,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Note: This item requires accessories for use | [Learn More](#)

Get Product Downloads

Diameter (mm):
160.00

Length (mm):
329.50

Type of Illumination:
Specialty

Note:
Aperture is not adjustable

Manufacturer: Edmund Optics®
Geometry: Backlight
Beam Diameter (mm): 138
Mounting Flange Diameter (mm): 200
Light Mount Inner Diameter (mm): 8

Regulatory Compliance

RoHS: Compliant
Certificate of Conformance: View
