

TECHSPEC® 182mm Telecentric Backlight Illuminator



182mm Telecentric Backlight Illuminator, #35-426

Stock **#35-426** 3-4 DAYS



1



€3.235,00

ADD TO CART

Qty 1+

€3.235,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Note: This item requires accessories for use | [Learn More](#)

Get Product Downloads

Diameter (mm):
208.00

Length (mm):
385.00

Type of Illumination:
Specialty

Note:
Aperture is not adjustable

Manufacturer: Edmund Optics®
Geometry: Backlight
Beam Diameter (mm): 182
Mounting Flange Diameter (mm): 260
Light Mount Inner Diameter (mm): 8

Regulatory Compliance

RoHS: Compliant
Certificate of Conformance: View
