

20 x 20mm Opal Checkerboard Target



Opal Checkerboard Calibration Target

Stock #12-197 3-4 DAYS



1



€482,00

ADD TO CART

Qty 1-4

€482,00

Qty 5+

€458,60

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Coating:
Blue Chrome (reflectivity <7% @ 550nm)

Substrate:
Opal Glass

Pattern Size (mm):
0.2 x 0.2 ±0.002

Dimensions (mm): 20 x20 ±0.10
Thickness (mm): 1.52 ±0.10
NIST Certification: No

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 235: Compliant
