

TECHSPEC®

242mm Telecentric Backlight Illuminator



TECHSPEC® Telecentric Backlight Illuminators

Stock #12-229 3-4 DAYS



1



€5.150,00

ADD TO CART

Qty 1+

€5.150,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Diameter (mm):
271.00

Length (mm):
678.00

Type of Illumination:
Specialty

Note:
Aperture is not adjustable

Manufacturer:

Edmund Optics®	
Geometry:	Backlight
Beam Diameter (mm):	242
Mounting Flange Diameter (mm):	298.5
Light Mount Inner Diameter (mm):	8

Regulatory Compliance

RoHS:	Compliant
Certificate of Conformance:	View
