


In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT



25 - 30° C Temp Range, Liquid Crystal Dot Sticker (30/Pack)

Stock #16-782 3-4 DAYS



1



€16,80

ADD TO CART

Qty 1-5

€16,80

Qty 6+

€14,70

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Operating Temperature (°C):
+25 to +30

Operating Temperature (°F):
77 - 86

Diameter (inches):
1.00

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
