

In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union



SUBMIT


TECHSPEC® 2mm Aperture, 75mm Length, Light Pipe Homogenizing Rod Mount



6mm Aperture, 50mm Length, Light Pipe Homogenizing Rod Mount, #64-913



Stock #64-906 3-4 DAYS



€115,00

ADD TO CART

Qty 1-5

€115,00

Qty 6+

€92,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Length of Compatible Light Pipe (mm):
75.00

Aperture of Compatible Light Pipe (mm):
2.00

Regulatory Compliance

RoHS:
Not Compliant

Certificate of Conformance:
[View](#)