

TECHSPEC® 400 - 750nm, 25mm Dia., Dichroic Longpass Filter Kit (8 Filters)



Stock **#15-102** [CONTACT US](#)

- 1 + **€1.215,00**

ADD TO CART

Qty 1+
€1.215,00

Volume Pricing
Request Quote

! Prices shown are exclusive of VAT/local taxes

Get Product Downloads



Diameter (mm): 25.00
Filters Included : Cut-on Wavelengths: 400nm , 450nm , 500nm , 550nm , 600nm , 650nm , 700nm , 750nm
Type: Longpass Filter
Number of Filters: 8

Environmental Durability:
Humidity per ML-STD-810H, Section 507.6
Temperature per ML-STD-810H, Section 501.7 and 502.7

Physical Durability:
Adhesion per ML-PRF-13830B, Section C.4.5.12
Moderate abrasion per ML-PRF-13830B, Section C.4.5.11
Cleaning per ML-C-48497A Section 4.5.4.2

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
