

50 x 50mm Opal Checkerboard Target



Opal Checkerboard Calibration Target

Stock #12-200 3-4 DAYS



1



€654,00

ADD TO CART

Qty 1-4

€654,00

Qty 5+

€621,40

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Coating:
Blue Chrome (reflectivity <7% @ 550nm)

Substrate:
Opal Glass

Pattern Size (mm):
0.8 x 0.8 ±0.002

Dimensions (mm): 50 x 50 ±0.10
Thickness (mm): 1.52 ±0.10
NIST Certification: No

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 235: Compliant
