




TECHSPEC Multi-Element Outer Tubes

Stock #37-902 3-4 DAYS

$\ominus$ 1 $\oplus$

€21,00

ADD TO CART

Qty 1-4	
€21,00	

Qty 5+
€18,70

Volume Pricing

[Request Quote](#)

! Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Function: Extension Tube (Fixed)
Thread Type: M29 x 0.75
Extension Length (mm): 50.0

Total Length (mm):
50.0

Construction:
Matte Black Anodized Aluminum (6061-T6)

Outer Diameter (mm):
32.0

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)