


In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

TECHSPEC® 80mm Length, Outer Tube



TECHSPEC Multi-Element Outer Tubes

Stock #37-903 3-4 DAYS



1



€28,40

ADD TO CART

Qty 1-4

€28,40

Qty 5+

€25,30

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Function:
Extension Tube (Fixed)

Thread Type:
M29 x0.75

Extension Length (mm):
80.0

Total Length (mm):
80.0

Construction:
Matte Black Anodized Aluminum (6061-T6)

Outer Diameter (mm):
32.0

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)