

In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union



SUBMIT

TECHSPEC® Liquid Lens Holder for 50mm Cx Lens



#33-686: Liquid Lens Holder for 50mm Cx Lens

Stock #33-686 3-4 DAYS



1



€25,00

ADD TO CART

Qty 1+

€25,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Lens Accessory

Compatible Lens:
[#33-566](#)

Lens Wavelength Range:
VIS

Regulatory Compliance

RoHS:

Compliant

Certificate of Conformance:
[View](#)