

In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union



SUBMIT

TECHSPEC® Multi-Element Retaining Ring Set



Multi-Element Retaining Ring Set

Stock #37-774 3-4 DAYS



1



€17,90

ADD TO CART

Qty 1-4

€17,90

Qty 5+

€15,70

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Function:

Retainer Ring

Thread Type:

M29 x0.75

Construction:

Matte Black Anodized Aluminum (6061-T6)

Thickness (mm): 3.00
Inner Diameter (mm): 25.0

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant

