

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT



Norland Optical Adhesive NOA73, 100g Bottle

Stock #17-345 3-4 DAYS



1



€119,00

ADD TO CART

Qty 1-4

€119,00

Qty 5+

€113,30

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Size (oz):
3.52

Glass Bonding:
Excellent

Metal Bonding:
Good

Plastic Bonding: Fair
Viscosity (cps): 130
Durability: Flexible
Index of Refraction n_d: 1.56 @ 589nm
Norland Number: 73
Shelf Life: 4 months
Type: Bottle
Typical Applications: Bonding delicate parts or parts where stress may be encountered
Absorption Range (nm): 350 - 380
Bonding Type: Glass to Glass
Energy for Full Cure (J/cm²): 4

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 240: Compliant
