

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

Pedestal Extender 11.0mm, M4, 1/4-20



Pedestal Extender

Stock #16-586 3-4 DAYS

- 1 + €9,90

ADD TO CART

Qty 1-4

€9,90

Qty 5+

€9,00

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
M4x0.7, 1/4-20

Type:
Metric/English

Construction:
Stainless Steel

Length (mm): 11.00
Function: Pedestal Extender
Post Diameter (mm): 25.4

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant