

Pedestal Extender 11.0mm, M4, M6



Pedestal Extender

Stock #16-312 3-4 DAYS

-

1

+

 €9,90

ADD TO CART

Qty 1-4
€9,90

Qty 5+
€9,00

Volume Pricing
Request Quote

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
M4x0.7, M6x1.0

Type:
Metric

Construction:
Stainless Steel

Length (mm): 11.00
Function: Pedestal Extender
Post Diameter (mm): 25.4

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant