

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

Pedestal Extender 20.8mm, M6, M4



Pedestal Extender

Stock #16-459 3-4 DAYS



1



€10,80

ADD TO CART

Qty 1-4

€10,80

Qty 5+

€9,70

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
M4x0.7, M6x1.0

Type:
Metric

Construction:
Stainless Steel

Length (mm): 20.80
Function: Pedestal Extender
Post Diameter (mm): 25.4

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant