

Pedestal Post, English, 12.7mm



Pedestal Post

Stock #15-845 3-4 DAYS

⊖ 1 ⊕ **€22,00**

ADD TO CART

Qty 1-4
€22,00
Qty 5+
€19,90
Volume Pricing
Request Quote

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
8-32

Type:
English

Weight (kg):
0.06

Construction: Stainless Steel
Length (mm): 12.70
Note: With 1/4-20 screw holes on the base
Function: Pedestal Post
Post Diameter (mm): 25.4
Base Diameter (mm): 31.5

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant