

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

Pedestal Post, Metric, 50.8mm



Pedestal Post

Stock #15-842 [CONTACT US](#)



1



€26,40

ADD TO CART

Qty 1-4

€26,40

Qty 5+

€23,80

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
M4x0.7

Type:
Metric

Weight (kg):
0.19

Construction: Stainless Steel
Length (mm): 50.80
Note: With M6x1.0 screw holes on the base
Function: Pedestal Post
Post Diameter (mm): 25.4
Base Diameter (mm): 31.5

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant
