

Pedestal Post, Metric, 76.2mm



Pedestal Post

Stock **#15-843** 3-4 DAYS

⊖ 1 ⊕ **€30,80**

ADD TO CART

Qty 1-4
€30,80

Qty 5+
€27,80

Volume Pricing
Request Quote

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
M4x0.7

Type:
Metric

Weight (kg):
0.31

Construction: Stainless Steel
Length (mm): 76.20
Note: With M6x1.0 screw holes on the base
Function: Pedestal Post
Post Diameter (mm): 25.4
Base Diameter (mm): 31.5

Regulatory Compliance

RoHS 2015: Compliant
Certificate of Conformance: View
Reach 233: Compliant
