


In order to provide more accurate inventory information, please select your shipping country:

Select Your Country/Region: European Union



SUBMIT

Resin 1.2mm Tip Width, Non-Marring Tweezers



Resin 1.2mm Tip Width, Non-Marring Tweezers, #36-449

Stock #36-449 3-4 DAYS



1



€30,00

ADD TO CART

Qty 1-9

€30,00

Qty 10+

€26,50

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Length (mm):
125.00

Type:
Resin Tweezer

Weight (g):
8

Tip Width (mm): 1.2
Tip to Tip Distance, Open (mm): 10
Holding Strength (g): 60

Regulatory Compliance

RoHS: Compliant
Certificate of Conformance: View
