


10mm Scale/100Div, 21mm Dia. Microscope Eyepiece Reticle



Stock #30-075

CLEARANCE

12 In Stock

-

1

+

€251⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-4	€251,00 each
Qty 5-9	€239,30 each
Qty 10+	€225,20 each
Need More?	Request Quote

!

 Prices shown are exclusive of VAT/local taxes

Product Downloads

Physical & Mechanical Properties	
21.00 +0.0/-0.1	Diameter (mm):
1.50 ±0.1	Thickness (mm):
±3	Overall Accuracy (µm):
	Centering (mm):

0.1	
	Line Width (mm):
0.012	
	Line Width Tolerance (mm):
±0.002	
Regulatory Compliance	
	RoHS 2015:
Compliant	
	Certificate of Conformance:
View	
	Reach 240:
Compliant	

Product Details

- Designed for 1XObjective Power
- 5 or 10mm Scale
- 1.5mm Thickness

Crossed Mcrometer Scales Transmission Reticles are clear glass discs marked with scales or patterns, which fit into microscope eyepieces at the internal real image plane. This is the intermediate image of the specimen/object under inspection produced by the objective lens of the microscope. By placing the reticle at this position, the scale/pattern is superimposed on the image, and both are in simultaneous focus. Crossed Mcrometer Scales Transmission Reticles are all designed for an objective power of 1X. These reticles are available in a 19 or 21mm diameter.