

**TECHSPEC® 12.7mm - 12mm (Fixed), Right Angle Post Clamp**



#53-357



Stock #58-671 **20+ In Stock**

- 1 + €28<sup>00</sup>

**ADD TO CART**

Volume Pricing

Qty 1-9	€28,50 each
Qty 10+	€25,70 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

**General**

Universal

Type:

Fixed

Function:

**Physical & Mechanical Properties**

