

[See all 60 Products in Family](#)

TECHSPEC® 30mm, Center Drive, Thru-Hole, 0.25" Travel, Metric Micrometer



TECHSPEC® Crossed-Roller Bearing Linear Translation Stages (Thru-Hole Model)



Stock #83-509 **5 In Stock**

⊖ 1 ⊕ €475.⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-2	€475,00 each
Qty 3+	€414,90 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Metric **Type:**

Physical & Mechanical Properties

Linear (X) **Type of Movement:**

