

Type-A to Micro D-Sub , USB 2.0 Locking Cable with Trigger, 5m



#59-390

Stock **#59-390** 3-4 DAYS



1



€109,00

ADD TO CART

Qty 1+

€109,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Camera Accessory

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)

Reach 233:

