

TECHSPEC® 70mm, Side Drive, Thru-Hole, English Extended Micrometer



Stock #66-443 **2 In Stock**

⊖ 1 ⊕ €640⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-2	€640,00 each
Qty 3+	€558,20 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General

English **Type:**

Physical & Mechanical Properties

Linear (X) **Type of Movement:**

Crossed Roller Bearing **Guide System:**

