

TECHSPEC® 90° Kinematic Adapter



Stock #34-248 3-4 DAYS

⊖ 1 ⊕

€129,00

ADD TO CART

Qty 1+
€129,00

Volume Pricing
[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Function:
Angle Adjustment

Length (mm):
62.40

Width (mm):
76.40

Thread Type:
C-Mount

Bolt Circle Diameter for M3 Tapped (mm): 38.0
Compatible Post: M6 x 1.0
Note: 25.0/25.4mm TECHSPEC® Kinematic Mount required
Construction: Black Anodized Aluminum

Regulatory Compliance

RoHS 2015: Compliant
Reach 224: Compliant
Certificate of Conformance: View
