

TECHSPEC®

Multi-Element Tube Spanner Wrench



Multi-Element Tube Spanner Wrench

Stock #11-046 3-4 DAYS

⊖

1

⊕

€52,50

ADD TO CART

Qty 1+

€52,50

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Function:
Spanner Wrench

Total Length (mm):
120.0

Outer Diameter (mm):
28.0

Inner Diameter (mm):
15.0

Regulatory Compliance

RoHS 2015:
[Compliant](#)

Certificate of Conformance:
[View](#)
