

In order to provide more accurate inventory information, please select your shipping country.


Select Your Country/Region: European Union



SUBMIT

Olympus Wide Field Tube Lens to T-Mount Holder

See More by [Olympus](#)



Olympus Wide Field Tube Lens to T-Mount Holder (#11-149)

Stock #11-149 3-4 DAYS



1



€149,00

ADD TO CART

Qty 1+

€149,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Mounting Threads:
External: M42 x0.75 (Male and Female)
Internal: M41 x0.5

Outer Diameter (mm):
56.0

Note:
Requires adapter [#52-301](#) to connect Olympus objective

Manufacturer:

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
