


In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

Replacement Tips for Resin 1.2mm Tip



Resin 1.2mm Tip Width, Non-Marring Tweezers, #36-449



Stock #36-450 3-4 DAYS



1



€15,00

ADD TO CART

Qty 1+

€15,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Length (mm):
48.00

Type:
Resin Tip Replacement

Tip Width (mm):
1.2

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
