

TECHSPEC[®]Aperture Kit for Achromat Prototyping Kit



Aperture Kit for Achromat Prototyping Kit

Stock **#54-720** **1 In Stock**

-

1

+

€156⁵⁶

ADD TO CART

Volume Pricing	
Qty 1-5	€156,56 each
Qty 6-10	€145,64 each
Qty 11-25	€138,84 each
Need More?	Request Quote

! Prices shown are exclusive of VAT/local taxes

Product Downloads

SPECIFICATIONS

General

Prototyping Kit

Type:

