

TECHSPEC[®] Cell with Two Retainers and Apertures for Achromat Proto Kit



Stock **#54-719** [CONTACT US](#)

-1+

€300^{.00}

ADD TO CART

Volume Pricing	
Qty 1-5	€300,00 each
Qty 6-10	€264,60 each
Qty 11-25	€244,00 each
Need More?	Request Quote

Prices shown are exclusive of VAT/local taxes

Product Downloads

General	
Prototyping Kit	Type:
Kit	Function:

Regulatory Compliance

