

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

TECHSPEC® Cell with Two Retainers for Achromat Prototyping Kit



Cell with Two Retainers for Achromat Prototyping Kit

Stock #58-255 ■ In Stock



1



€146,10

ADD TO CART

Qty 1-5

€146,10

Qty 6+

€128,70

Volume Pricing

[Request Quote](#)

ⓘ Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Type:
Prototyping Kit

Function:
Kit

Regulatory Compliance

RoHS:
[Compliant](#)

Certificate of Conformance:
[View](#)
