

# Clear, 10cc Barrel (100 count), Dispenser Barrel



Clear Dispenser Barrel

Stock **#25-096** 1 In Stock

-

1

+

€279<sup>00</sup>

ADD TO CART

| Volume Pricing |                               |
|----------------|-------------------------------|
| Qty 1+         | €279,00 each                  |
| Need More?     | <a href="#">Request Quote</a> |

Prices shown are exclusive of VAT/local taxes

## Product Downloads

Regulatory Compliance

Certificate of Conformance:

[View](#)

## Product Details

- Empty Barrels are Available in Clear, Amber, and Opaque Versions
- Tips are Available in Plastic and Steel

- [Dispensing Guns](#) and [Pre-Loaded Norland Adhesive Barrels](#) are Also Available

Dispenser barrels, tips, and other accessories are available for use with dispensing guns and hand plungers. Dispenser Barrels are available in a clear version for normal fluids, an amber version for UV-curing fluids, and an opaque version for light sensitive fluids. Dispenser tips are available in plastic and steel to ensure varying levels of controllable deposits without leakage.

**Note:** Barrels, tips, and adhesive are sold separately from [dispensers](#). [Pre-loaded Norland adhesive barrels](#) are also available.