

TECHSPEC[®] Complete Achromat Prototyping Kit



Stock **#54-666** [CONTACT US](#)

-

1

+

€1.270⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-5	€1.270,00 each
Qty 6-10	€1.118,60 each
Qty 11-25	€1.029,60 each
Need More?	Request Quote

!

 Prices shown are exclusive of VAT/local taxes

Product Downloads

General	
Prototyping Kit	Type:
Kit	Function:

Regulatory Compliance

