

[See all 70 Products in Family](#)

**TECHSPEC® DAPI (High Contrast) Filter Cube Set - Nikon**



Olympus Filter Cube #86-833



Stock #86-371-NKN [CONTACT US](#)

⊖ 1 ⊕ €1.350<sup>00</sup>

**ADD TO CART**

Volume Pricing	
Qty 1-9	€1.350,00 each
Qty 10+	€1.215,00 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

**General**

Fluorescence Filter Kit **Type:**

DAPI (High Contrast) **Compatible Fluorophore:**

[#84-094](#) **Excitation Filter:**

[#84-095](#) Emission Filter:

[#86-330](#) Dichroic Filter:

Compatible Microscope Brand:  
Nikon

## Optical Properties

[Fused Silica](#) (Coming 7980) Substrate:

Hard Coated Coating:

Excitation Wavelength (nm):  
382 - 393

Emission Wavelength (nm):  
417 - 477

Dichroic Cut-On Wavelength (nm):  
409.00

Wavelength Range (nm):  
382 - 477

## Regulatory Compliance

Certificate of Conformance:  
[View](#)

## Product Details

- Ensure Proper Alignment of Filter Sets
- Easily Adapt to Common Microscopes
- Available in a Variety of Wavelengths
- Unmounted [TECHSPEC® Fluorescence Filter Kits](#) Also Available

TECHSPEC® Pre-Mounted Fluorescence Filter Cube Sets are designed to easily integrate into popular Olympus, Nikon, Zeiss, and Zaber microscopes. Each fluorescence filter cube set consists of one [TECHSPEC Fluorescence Dichroic Filter](#) and two [TECHSPEC Fluorescence Bandpass Filters](#) that have been specifically oriented to achieve optimum excitation and emission of popular fluorophores. Additionally these filters are also used in polymerase chain reaction (PCR) diagnostic instruments. For increased adaptability, filters can be easily removed and replaced to meet the needs of the application. Unmounted [TECHSPEC Fluorescence Filter Kits](#) are also available.

**Built-to-Order Product:** Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).

**Built-to-Order Product:** Each product is assembled from stock components at the time of your order. If all components are in-stock, your product will typically ship within 3 - 5 days. If you have any questions about delivery or inventory status, please [contact us](#).