

[See all 6 Products in Family](#)

## Diffuse reflectance phantom Adipose Tissue 50 x 50 x 10 mm



Stock #71-280 **2 In Stock**

⊖ 1 ⊕ €1.175<sup>00</sup>

**ADD TO CART**

Volume Pricing	
Qty 1-4	€1.175,00 each
Qty 5-9	€1.057,50 each
Need More?	<a href="#">Request Quote</a>

**i** Prices shown are exclusive of VAT/local taxes

### Product Downloads

### General

Adipose Tissue **Type:**

### Physical & Mechanical Properties

50 x 50 x 10 **Dimensions (mm):**

### Material Properties

**Reduced Scattering Coefficient,  $\mu'_s$ :**

Absorption Coefficient,  $\mu_a$ :

0.1

---

## Regulatory Compliance

Certificate of Conformance:

[View](#)

---

## Product Details

- Mimic Muscle, Brain, or Adipose Tissue
- Available in 25 x 25 mm and 50 x 50 mm Sizes
- Ideal for a Range of Life Science Applications

Diffuse Reflectance Phantoms are solid samples that mimic the optical properties of biological tissues for testing, calibration, and characterization of biophotonics devices. These homogeneous phantoms, with well characterized optical properties at 780nm, simulate the optical properties and mechanical characteristics of biological tissues such as muscle, brain, and adipose tissue. These samples come in two sizes, 25 x 25mm or 50 x 50mm, to accommodate multiple experimental setups. Diffuse Reflectance Phantoms are ideal for spatial frequency domain imaging, diffuse reflectance spectroscopy, diffuse Raman spectroscopy, and biophotonics/ biomedical imaging.

---