

Diffuse reflectance phantom Brain 50 x 50 x 10 mm



Stock **#71-278** 2 In Stock

-

1

+

€1.175⁰⁰

ADD TO CART

Volume Pricing	
Qty 1-4	€1.175,00 each
Qty 5-9	€1.057,50 each
Need More?	Request Quote

Prices shown are exclusive of VAT/local taxes

Product Downloads

SPECIFICATIONS

General

Type:

Brain

Physical & Mechanical Properties

50 x 50 x 10	Dimensions (mm):
Material Properties	
12	Reduced Scattering Coefficient, μ_s' :
0.2	Absorption Coefficient, μ_a :
Regulatory Compliance	
View	Certificate of Conformance:

PRODUCT DETAILS

- Mimic Muscle, Brain, or Adipose Tissue
- Available in 25 x 25 mm and 50 x 50 mm Sizes
- Ideal for a Range of Life Science Applications

Diffuse Reflectance Phantoms are solid samples that mimic the optical properties of biological tissues for testing, calibration, and characterization of biophotonics devices. These homogeneous phantoms, with well characterized optical properties at 780nm, simulate the optical properties and mechanical characteristics of biological tissues such as muscle, brain, and adipose tissue. These samples come in two sizes, 25 x 25mm or 50 x 50mm, to accommodate multiple experimental setups. Diffuse Reflectance Phantoms are ideal for spatial frequency domain imaging, diffuse reflectance spectroscopy, diffuse Raman spectroscopy, and biophotonics/biomedical imaging.