


[See all 5 Products in Family](#)

Lumedica OQ LabScope 3.0 /R OCT Imaging System



OCT Imaging System (OCT Stand, Monitor, Mouse, and Keyboard Sold Separately as #37-140)



Stock **#29-144** [CONTACT US](#)

-

1

+

€17.495⁰⁰

ADD TO CART

Volume Pricing	
Qty 1+	€17.495,00 each
Need More?	Request Quote

! Prices shown are exclusive of VAT/local taxes

Note: This item requires accessories for use | [Learn More](#)

Product Downloads

General

Model Number:

OQ LabScope 3.0/R

A-Scan Line Rate (lines/second):

34000

B-Scan Image Rate (lines/second) :

Physical & Mechanical Properties	
6	Weight (lbs):
2.72	Weight (kg):
Dimensions (mm): System: 413 x216 x 153 Scanner: 41 x 172 x 67	
Depth Resolution (µm): Air: 5 Tissue: 3	
Transverse Resolution (µm): 18	
Volume Scan Range (mm, diameter of a circle): 5mm	
Optical Properties	
Center Wavelength CWL (nm): 860.00	
Linear Scan Range (mm): 7mm	
Axial Resolution (In Air / In Tissue): 5 µm / 3 µm	
Axial Scan Range (Maximum Imaging Depth in Air / In Tissue): 2mm / 1.3mm	
Sensor	
Image Size (Pixels): 512 x 512	
Electrical	
Sensitivity (OSNR): 100 dB	
Output Power (µW): 750	
Regulatory Compliance	
Certificate of Conformance: View	

Product Details

- Affordable Optical Coherence Tomography Imaging System
- Ideal for Biological Sample Imaging, Sample Characterization, and OCT Education
- Compact Form Factor Tabletop Device

The Lumedica OCT Imaging System is an all-in-one high resolution Optical Coherence Tomography Imaging System. Designed as an affordable alternative to prohibitively expensive and complex units, the system provides competitive resolution, frame-rate, and power specifications at an industry-leading price. With a laser centered at 840nm, this OCT Imaging System is ideal for inspection at depths of up to 2.8mm in air, and to 2mm in tissue. The Lumedica OCT Imaging System features axial resolutions of 5µm and lateral resolutions of 15µm, matching the performance of much more expensive units. With a compact form factor, integrated computer, and intuitive graphical user interface for laboratory or field use, the Lumedica OCT Imaging System is an ideal tool for users of all experience levels. An optional workstation package includes a scanner stand with vertical stage, large monitor, as well as a wireless mouse/keyboard combination to interface with the device.