

[See all 5 Products in Family](#)

M21 x 0.5 Eyepiece Retainer Ring for #35-689



M21 x 0.5 Eyepiece Retainer Ring for #35-689

Stock **#17-409** **20+ In Stock**

⊖ 1 ⊕ €23.⁰⁰

ADD TO CART

Volume Pricing	
Qty 1+	€23,00 each
Need More?	Request Quote

ⓘ Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Accessory **Type:**

Compatible with [#35-689](#) **Note:**

Physical & Mechanical Properties

17 **Clear Aperture CA (mm):**

Thickness (mm):

3.00 ±0.05

Black Oxidized Steel **Construction:**

Threading & Mounting

M21.0 x 0.5 **Thread Type:**

Regulatory Compliance

Compliant **RoHS 2015:**

View **Certificate of Conformance:**

Compliant **Reach 235:**


Product Details

- Widefield and DIN-45mm Versions Available
- Widefield Versions Offer Adjustable Diopter
- Three Standard Magnifications

Microscope Eyepieces attach to a variety of microscope systems featuring designs that consist of an achromatic doublet with a single plano-convex lens, where both convex surfaces face each other. These eyepieces are ideal for use with standard achromatic objectives and upright, brightfield microscope systems. Microscope Eyepieces are available as standard parts with different magnifications, focal lengths, and field stops. Additional magnification options available upon request - [contact us](#) today to discuss your specific requirements.

Note: Eyepieces are sold individually, not as a set.

Technical Information



Stock No.	A	B	C	D
#35-689	30.1	24.5	10	28
#35-690	33.7	26.3	10	42
#35-691	19.6	26.3	10	40
#35-692	19.5	24.5	10	28