

[See all 12 Products in Family](#)

**TECHSPEC® Mounting Clamp, 2.25" Centerline, 30mm ID**



Mounting Clamp, 2.25" Centerline, #56-871



Stock **#56-871** **20+ In Stock**

− 1 + €190<sup>00</sup>

**ADD TO CART**

Volume Pricing

Qty 1+	€190,00 each
Need More?	<a href="#">Request Quote</a>

Prices shown are exclusive of VAT/local taxes

Product Downloads

General

Stock No. of Mounting Clamp:  
[#56-870](#) Sold Separately

Accessory **Type:**

Physical & Mechanical Properties

Inner Diameter (mm):

