

TECHSPEC® Mounting Clamp for 2 - 8X Research-Grade Variable Beam Expander



Mounting Clamp for 2 - 8X Research-Grade Variable Beam Expander, #89-295

Stock #89-295 3-4 DAYS



1



€164,00

ADD TO CART

Qty 1-4

€164,00

Qty 5+

€145,00

Volume Pricing

[Request Quote](#)

⚠ Prices shown are exclusive of VAT/local taxes

Get Product Downloads



Type:
Mounting Clamp

Regulatory Compliance

RoHS:
[Compliant](#)

