

# Olympus MPLN 20X Objective

See More by [Olympus](#)



Stock #87-137 **NEW** 2 In Stock

⊖ 1 ⊕ €1.150<sup>00</sup>

**ADD TO CART**

### Volume Pricing

Qty 1+	€1.150,00 each
Need More?	<a href="#">Request Quote</a>

ⓘ Prices shown are exclusive of VAT/local taxes

### Product Downloads

### General

Model Number:  
MPLN20X

Compatible Tube Lens Focal Length (mm):  
Focal Length: 180mm

Type:  
Microscope Objective

Style:  
Infinity Corrected

Olympus

Manufacturer:

## Physical & Mechanical Properties

1.1 Field of View (mm):

43.70 Length excluding Threads (mm):

24 Maximum Diameter (mm):

111 Weight (g):

## Optical Properties

N/A Compatible Cover Glass Thickness (mm):

9.00 Focal Length FL (mm):

20X Magnification:

0.40 Numerical Aperture NA:

0.84 Resolving Power ( $\mu\text{m}$ ):

1.72 Depth of Field ( $\mu\text{m}$ ):

1.3 Working Distance (mm):

400 - 700 Wavelength Range (nm):

22 Field Number (mm):

45 Parfocal Length (mm):

N/A Immersion Liquid:

8.10 Entrance Pupil Diameter (mm):

## Threading & Mounting

RMS / 20.32mm x 36 TPI Mounting Threads:

## Regulatory Compliance

[View](#) Certificate of Conformance:

## Product Details

- Ideal for Transmitted or Reflected Brightfield Microscopy
- Produce Flat Images Up to Field Number (F.N.) 22
- Excellent Price-to-Performance Ratio
- Additional [Olympus Infinity Corrected Objectives](#) Available

Olympus MPlan Achromatic Objectives are designed for brightfield observations, providing excellent contrast and optimum flatness throughout the field of view. Engineered for inspection of silicon wafers, metals, and other industrial microscopy applications, these objectives may also be used in Raman microscopes and other reflected light techniques without coverslip correction. Olympus MPlan Achromatic Objectives are available in 5X-100X magnifications and provide excellent image flatness up to F.N. 22. [Olympus Plan Achromatic Objectives](#) are also available for observing biological specimens.