

In order to provide more accurate inventory information, please select your shipping country.

Select Your Country/Region: European Union



SUBMIT

TECHSPEC® Polarization Tester



Polarization Tester

Stock #37-699 3-4 DAYS



1



€36,00

ADD TO CART

Qty 1+

€36,00

Volume Pricing

[Request Quote](#)

Prices shown are exclusive of VAT/local taxes

Get Product Downloads

Length (mm):  
200.00

Thickness (mm):  
0.75

Transmission (%):  
42

Type:  
Polarization Tester

Construction:

|                                            |
|--------------------------------------------|
| Polarizing Film                            |
| <b>Wavelength Range (nm):</b><br>400 - 700 |
| <b>Width (mm):</b><br>20.00                |

Regulatory Compliance

|                                                            |
|------------------------------------------------------------|
| <b>RoHS 2015:</b><br><a href="#">Compliant</a>             |
| <b>Reach 224:</b><br><a href="#">Compliant</a>             |
| <b>Certificate of Conformance:</b><br><a href="#">View</a> |

---