

[See all 12 Products in Family](#)

## T-Thread Mounting Ring for Motorized Filter Wheel



#59-772: T-Thread Mounting Ring for Motorized Filter Wheel

Stock **#59-772** **2 In Stock**

1  €107<sup>00</sup>

**ADD TO CART**

### Volume Pricing

Qty 1+	€107,00 each
Need More?	<a href="#">Request Quote</a>

**!** Prices shown are exclusive of VAT/local taxes

### Product Downloads

### Regulatory Compliance

[Compliant](#)

RoHS 2015:

[Compliant](#)


Reach 224:

[View](#)

Certificate of Conformance:



### Top View of Mounting Plates



### Side View of Mounting Plates

